



Ultrasonic, vibrating facial skin scrubber

Generuota 2025-04-03

Name

Description

More

Ultrasonic Vibrating FacialExfoliator Main functions: * Deep cleansing: helps deeply cleanse the skin through high frequency vibration to remove blackheads, acne, dead skin, make-up residue and pore oils. High frequency vibration breaks down melanin and dark spots, leaving your skin smooth, firm and supple. * Improved absorption of skincare essence: negative ions affect cellular activity, a facial scrub can increase your skin's ability to absorb nutrients and cream.

Material: ABS+304 stainless steel

Battery: 3. 7V/550mah

Vibration frequency: 2500Hz/S

Output power: 2w

Product size: 195*76*38mm

Weight: 86g

Charging time: about 3hrs

Usage time after full replacement: about 3hrs

3 models: cleansing/ionic input/EMS lifting mode

Packing list: 1 pcs Ultrasonic Skin Scrubber, 1 pcs USB cable, 1 pcs user manual.

FUNCTIONS

Step 1 - Cleansing. Ultrasonic cleansing mode - PEELING (for use with facial cleanser/facial cleansing milk). Wet your face with warm water and apply a facial cleanser/facial cleansing milk that foams. Press the button on the switch on the device, then select the PEELING mode. With the spatula pointing inwards and the head of the device downwards, press the spatula against the skin and rub the head into the pores for 3 minutes.

Step 2 - Deep cleansing. Positive ion output mode - ION+ (for use with toner/lotion). Apply a toner-lotion to the face. Press the button briefly to select ION+ mode. With the head pointing inwards, press the head against the skin, then lower the head towards the pores to deep cleanse, do this for 3 minutes. Then rinse your face with warm water.

Step 3 - Nutrient infusion. Negative ion output mode - ION- (use with essence/lotion/mask). Apply a suitable serum/lotion to the face. Select ION mode. Pulling out the spade with the hooves pointing upwards, bring it to the skin, then tilt the spade head towards the pores and massage for 3 minutes.

Step 4 - Soothing massage. EMS Pulse Mode-EMS PULSE (with firming cream). Select EMS PULSE mode, pull out the paddle, bring it close to the skin, then point the head of the device towards the pores and massage to soothe the skin for 3 minutes.



SMART ELECTRONICS
SALE / REPAIR



CONTACT US | +370 639 77277

QUICK PRODUCTS DELIVERY

1-2 w.d.

Name	Description
Galia	2W
Colour	White
DIMENSIONS	195 * 76 * 16 * 38mm.
Battery	550 mAh

UAB "REGO TECH"

Company code: 302449356
PVM code: LT100005005917
C/A: LT927300010118749098

WWW.MARSIIETIS.LT

Email: info@marsietis.lt
Mobile: +370 (639) 77 277
Vilnius st. 230 Siauliai, LT-76303, Lithuania