



### Wrinkle reduction vibrating machine with lights (black)

Generuota 2025-09-08

#### Name

#### Description

#### More

This product can be used to clean the skin, do EMS care and moisturize skin in a warm temperature, which integrates positive ion, negative ion, high-frequency vibration, red light, blue light and green light. Positive and negative ions: By utilizing the principle that opposite charges attract between positive and negative ions, derive and thoroughly clean the dirt from deep pores through bioelectric currents; The cosmetic components will penetrate into the deep skin through ion introduction. Vibration: Make human skin become tight and increase elasticity and increase the metabolism of the skin by adopting the methods of mechanical vibration and voltage stimulation. Tepidity: Tepid massage promotes absorption, stretches skin texture, and makes the skin become full.

How to use: Apply some skin care products on neck and face. Use the massage head close to neck, lift from bottom several times. Then fit the face skin and lift, make skin tightened. Lift along the wrinkles, remove double chin, smile lines, eye wrinkles and activates skin cell.

#### FUNCTIONS

Red light: it can effectively penetrate deep into the skin, revitalise collagen, remove darkness, reduce fine lines and wrinkles.

Blue light: it can facilitate collagen synthesis, rejuvenate the skin, shrink sagging skin, tighten skin pores and visibly improve oily, sensitive and acne-prone skin.

Green light: it can improve the cells' alternative oxygen utilisation function, facilitate skin microcirculation, accelerate lymphatic circulation and reduce oedema, and soothe the skin to a stable state.