



CONTACT US | +370 639 77277

QUICK PRODUCTS DELIVERY

1-2 w.d.



XPREEN Facial Pore Remover with LED display + Metal Pore Remover Tools

Generuota 2026-02-17

Name

More

Description

XPREEN facial pore cleaner with LED display. This cleanser is unique in that it cleans accumulated facial dirt and facial skin pores quickly and without any skin deformation or damage. removes acne, acne. This cleaner is designed based on a proven biomechanics principle, making it safe to use. The LED display allows you to select the desired mode and see in what mode the suction is working, which consists of three modes: weak mode, standard mode and powerful mode. The set consists of 4 different nozzles, which are adapted to different facial skin. It also provides a lot of convenience in using this device, as different nozzles allow access to different areas of the face such as: nose, eyebrow area, chin, areas around the ears. The use of this device not only results in a clean and glowing skin, but also a reduction in wrinkles, as the device has suction and stretching properties that can remove small facial wrinkles, making the face look younger. This machine is unique not only in its suction power and choice of different nozzles, but also in the set of metal pore removal tools that will allow you to care for and clean your skin at home.

NOTE: For best results when using this cleanser, it is recommended to warm the skin with steam or a warm cotton pad. The heat opens up the pores, making them easier to remove from the skin.

More information

Use twice a week. After the face is clean, smear the nose with a blackhead remover with a cotton pad and keep it for 3-5 minutes, you can use a hot air steam spray to fully open the pores. If there is no such possibility, a warm cotton pad will also work perfectly, and after holding it for a couple of minutes, the pores will open. Then you can already select the required head. After pressing the power button, hold it longer to turn on the device. First, choose a light cut and massage with light strokes to remove blackheads. Do not hold in the same place for more than 3 seconds. If necessary, the mode can be changed.



SMART ELECTRONICS
SALE / REPAIR



CONTACT US | +370 639 77277

QUICK PRODUCTS DELIVERY

1-2 w.d.

Name Description

ACCESSORIES 1 * cleaning machine

4 * tips

5 *metal pore cleaning tools

1 * USB cable

1 * user manual