

**CONTACT US | +370 639 77277** 

QUICK PRODUCTS DELIVERY

1-2 w.d.



## LED photonic light face mask, 7 colors

Generuota 2024-05-17

Name Description

Model LED photonic light face mask.

More

LED photonic light face mask, 7 colors. Let's first start with what phototherapy is: WHAT PHOTOTOTOTHERAPY IS - the technology of LED phototherapy application is the treatment of tissues with light of a selected wavelength (7 different options have been selected, for a description please see the auction). They penetrate the inside of the skin without direct intervention. The treatment is completely painless, stimulating various biological and chemical changes. The use of light greatly improves the condition of the skin and accelerates its regeneration. NON-INJURING AND EFFECTIVE TREATMENT - the use of the mask is in no way painful for our skin and its use has many benefits, namely the removal of acne and various discolourations or scars, the reduction of wrinkles, nourishment and improvement of the skin's elasticity, and delaying. The aging process also helps to get rid of puffiness under the eyes. Photon therapy using led mask is safe for all skin types. FLEXIBLE headband includes a flexible screen saver that allows you to adjust to the volume of your head. Phototherapy is recommended for both women and men, regardless of age. HOW DOES IT WORK? A beam of light penetrates the skin and reaches the cells inside, stimulating their processes. Its effectiveness can be compared to that of a laser, but the LED mask is painless to use. It does not affect the structure of the cells and naturally activates them. Remember! Before using the mask, wash off your make-up thoroughly. The mask illuminates only clean skin. For long-lasting effects, it is recommended to use 10-20 treatments every 2-3 days. They protect and preserve the eyes, no light passes through them, so the mask is completely safe and does not affect the eyes in any way. We have included an elastic band that will allow you to adjust to any head diameter!

## **FUNCTIONS**

Made of eco-friendly ABS material, non-toxic, odourless, healthy and durable. 7 red, blue, green, yellow, purple, white and cyan bulbs to choose according to your needs. Each bulb has 10 gears to adjust the intensity of the light. Each bulb has a timer function, with an effective control of 5/10/15/20/25/30 minutes. Stimulates skin metabolism, improves the appearance of blackheads, tightens pores. Cleanse the skin, break down pigment, dilute blemishes. Activate skin collagen secretion, dilute fine lines, firm skin. No surgery, no needle, safe and fast skin treatment and protection. The mask conforms to the contour of the face for a comfortable and wearable fit. Suitable for all skin types.

UAB "REGO TECH"

Company code: 302449356 PVM code: LT100005005917 C/A: LT927300010118749098 WWW.MARSIETIS.LT

Email: info@marsietis.lt Mobile: +370 (639) 77 277

Vilniaus st. 230 Siauliai, LT-76303, Lithuania



CONTACT US | +370 639 77277

QUICK PRODUCTS DELIVERY

1-2 w.d.

Name	Description
Colour	Apparatus: white. Number of colors/60/in: 7 + multicolor.
PAPILDOMA	Set the timer - 5/10/15/20/25/30/35/40/45/50/55min, choose the right colour - 7+ multicolours. You will see every setting on the transparent display! WS TREATMENT PRACTICES: acne, discolouration, wrinkles, dilated blood vessels, tired skin, irritation, small wounds, post-cleansing, needle and micro-needle mesotherapy.
Įtampa	220-240 V
VOLTAGE	Power: 10W / maximum power of 1 diode: 40mW.
ACCESSORIES	Note: Do not look directly at the LED light source for long periods of time. Do not read books, watch TV or surf the internet while using the photon mask. If you use a disposable mask or cosmetics with the LED mask, choose high quality products with natural ingredients. The mask material is medical grade silicon, which will smell a bit when you open the package, but the smell will disappear very quickly. Inside the silicon mask there are LED light sources and a PCB board. No folding, no squeezing, no strong pressure. Use the LED mask on the bed, and when you have finished the treatment, put the mask back on the box, do not leave it on the bed. Do not use water to clean the LED mask, do not drown in water. Wipe the LED mask with a damp cloth. Do not use any chemical liquid to clean the LED mask and soften the mask. Do not expose to direct sunlight Tips: The effect will be better if used with a high quality moisturising skin care product.

Company code: 302449356 PVM code: LT100005005917 C/A: LT927300010118749098 Email: info@marsietis.lt Mobile: +370 (639) 77 277

Vilniaus st. 230 Siauliai, LT-76303, Lithuania