



EMS foot electric smart pulse mat, massager

Generuota 2025-06-02

Name	Description
Model	EMS mat
More	The EMS Electric Impulse Mat / Massager is designed for those who wear high heels for long periods of time, have sore and tired feet, stand for long periods of time and their feet start to ache, freezing feet, old people with bleeding disorders, etc. The EMS is only able to massage the pads of the feet. Controlled by remote control, the palms of the hands can be supported. Relieve pain, fatigue and give you a good day!
More information	Product description : Name: EMS Size: 32 cm * 29 cm * 2,8 cm Power: USB Massage principle: vibration Mode intensity: 6 modes, 9 levels of intensity. Working voltage: 3.7V Working time: 15 minutes Working current: 15mA
Type	1. Place your feet flat on the massage mat, turn on the remote control and adjust it to your strength. 2. Press and hold M/+ for three seconds to switch the host computer on, switch off and press OFF /- to zero position, then press and hold M to switch off. 3. You can switch between six massage modes, which are: kneading, knocking, acupuncture, massage, acupressure, pulse. 4. Force adjustment. Press ON/+ to increase the intensity; press OFF/+ to decrease the intensity. 5. Recommended for 30 minutes each time, 1 - 2 times a day. 6. It is better to use the unit 30 minutes after a meal or bath.